

NASO-PHARYNGO-LARYNGOSCOPE

OP-20

HIGH RESOLUTION IN THE SMALLEST DIAMETER

- Pediatric use
- Great image quality
- High resistance



OPTOMIC®
producing quality

OPTOMIC flexible endoscopes provide excellent optical characteristics for the inspection of nasal passages, the pharynx and the larynx. This fiberscope has been designed to facilitate examination efficiency and enhance patient comfort.

NASO-PHARYNGO-LARYNGOSCOPE **OP-20**



FIBROLUX 250
LIGHT SOURCE



FIBROLUX LED
LIGHT SOURCE



FIBER OPTIC CABLE



FIBERSCOPE HOLDER




LEAKAGE TEST DEVICE



TWO-HANDED GRIP

	OP-20
Image guide system	11.000 pixels
Flexible part diameter	2.5mm
Depth of focus	5-50mm
Deflection angle	180° up / 90° down
Direction of view	0°
Field of view	90°
Working length	270mm

- High definition images
- Focusing adjustment.
- ETO valve.
- Soakable in disinfectant liquid.
- Light guide connector for Storz and ACMI (optional Wolf).
- Delivered in a protective carrying case.



producing quality

Specifications are subject to change without notice

