

# NASO-PHARYNGO-LARYNGOSCOPE

## OP-30

HIGH RESOLUTION IN THE SMALLEST DIAMETER

- Up to 18.000 pixels
- Detailed vision
- High resistance



OPTOMIC flexible endoscopes provide excellent optical characteristics for the inspection of nasal passages, the pharynx and the larynx. Both fiberscopes have been designed to facilitate examination efficiency and enhance patient comfort.

## NASO-PHARYNGO-LARYNGOSCOPE **OP-30**



FIBROLUX 250  
LIGHT SOURCE



FIBROLUX LED  
LIGHT SOURCE



FIBER OPTIC CABLE



FIBERSCOPE HOLDER




LEAKAGE TEST DEVICE



TWO-HANDED GRIP

	OP-30
Image guide system	18.000 pixels
Outer diameter	3,4mm
Flexible part diameter	3,2mm
Depth of field	2,5-50mm
Deflection angle	130°up/down
Direction of view	0°
Working length	300mm

- High definition images (up to 18.000 pixels).
- Focusing adjustment.
- ETO valve.
- Soakable in disinfectant liquid.
- Light guide connector for Storz and ACMI (optional Wolf).
- Delivered in a protective carrying case.



Specifications are subject to change without notice

